



### Product Description:

#### Epinephrine Injector Dispenser - 5 Injector Capacity

- Holds up to 5 epinephrine injectors
- Dual-sided clear PETG plastic door sleeve for identification signage
- Keyholes for wall mounting
- Made of quartz ABS and clear PETG plastic
- Made in U.S.A.

### Dispenses:

- Epinephrine auto injectors for anaphylactic emergencies (injectors not included).

### Product Specifications (Overall External Dimensions):

- 9"W x 11 <sup>3</sup>/<sub>16</sub>"H x 2 <sup>1</sup>/<sub>2</sub>"D (22.9cm x 28.4cm x 6.4cm)
- 1.85 lbs (0.8kg) - approximated

### Package Specifications:

- Package Qty 1 per case
- 14"W x 12"H x 9"D (35.6cm x 30.5cm x 22.9cm)
- 4.00 lbs (1.8kg) - approximated

### Mounting Specifications:

