



Product Description:

Glove Box Dispenser - Double - Space Saver

- Holds two boxes of gloves that fit end-to-end
- Two-way keyholes for vertical or horizontal wall mounting
- Made of white powder coated steel
- Made in U.S.A.

Dispenses:

- Most brand name standard glove boxes - Boxes no larger than 5 1/4" W x 10" H x 3 3/4" D (13.3cm x 25.4cm x 9.5cm)

Product Specifications (Overall External Dimensions):

- 5 5/8" W x 20 1/16" H x 3 15/16" D (14.2cm x 50.9cm x 10.0cm)
- 2.45 lbs (1.1kg) - approximated

Package Specifications:

- Package Qty: 4 per case
- 24" W x 16" H x 6" D (61.0cm x 40.6cm x 15.2cm)
- 14.00 lbs (6.4kg) - approximated

Mounting Specifications:

