



Product Description:

Glove Box Dispenser - Triple

- Holds three boxes of gloves
- Two-way keyholes for vertical or horizontal wall mounting
- Made of clear PETG plastic
- Made in U.S.A.

Dispenses:

- Most brand name standard glove boxes - Boxes no larger than 5 1/2"W x 10"H x 3 3/4"D (14.0cm x 25.4cm x 9.5cm)

Product Specifications (Overall External Dimensions):

- 15 15/16"W x 10 1/8"H x 3 15/16"D (40.5cm x 25.7cm x 10.0cm)
- 0.95 lb (0.4kg) - approximated

Package Specifications:

- Package Qty: 6 per case
- 21"W x 14"H x 10"D (53.3cm x 35.6cm x 25.4cm)
- 9.00 lbs (4.1kg) - approximated

Mounting Specifications:

