



Product Description:

Glove Dispenser

- Holds one box of gloves with a built-in spring
- Two-way keyholes for wall-mounting
- Made of 20 gauge stainless steel
- Made in USA

Dispenses:

Fits most brand name standard glove boxes (Boxes no larger than 6.44"W x 9.93"H x 3.74"D (16.4cm x 25.2cm x 9.5cm))

Product Specifications (Overall External Dimensions):

- 6.50"W x 9.96"H x 3.80"D (16.5cm x 25.3cm x 9.7cm)
- 2 lbs (.9 kg) - approximated

Package Specifications:

- Package Qty 12 per case
- 20"W x 16"H x 12"D (50.8cm x 40.7cm x 30.5cm)
- 24 lbs (11.9 kg) - approximated

Mounting Specifications:

