












How to Choose a Bowman® Protection Organizer

| Protection Organizer |  |  |  |  |  |  |  |  |  |
|--------------------------------|---|---|---|---|---|--|---|---|---|
| | LD-007 | LD-037 | LD-063 | LD-064 | LD-067 | LD-070 | PS007-0212 | PS012-0512 | PS015-0212 |
| Material | Quartz ABS Plastic | Quartz ABS Plastic | Quartz ABS Plastic | Quartz ABS Plastic | Quartz ABS Plastic | Quartz ABS Plastic | Quartz ABS Plastic | Quartz Powder Coated Aluminum | Quartz ABS Plastic |
| Dimensions | 16 5/8" W 25 1/4" H 4 1/2" D | 22 3/4" W 20 1/2" H 4 3/4" D | 19 11/16" W 30" H 4 3/4" D | 16 3/4" W 25 5/8" H 4" D | 19 11/16" W 28 5/8" H 4 3/4" D | 20 1/16" W 25 1/4" H 4 3/4" D | 20" W 20 15/16" H 4 15/16" D | 22 3/4" W 25 5/8" H 4 9/16" D | 25 1/8" W 30 3/4" H 4 3/4" D |
| Permanently Mounted / Portable | ● | ● | ● | ● | ● | ● | ● | ● | ● |
| 4" Deep Unit | | | | ● | | | | | |
| Reusable / Launderable Gowns | | ● | ● | ● | ● | | | ● | ● |
| Disposable Flat Pack Gowns | ● | ● | ● | ● | | ● | ● | ● | ● |
| Disposable Blue Gowns | | ● | ● | | | | ● | ● | ● |
| Gowns Single | ● | ● | | ● | | ● | | ● | |
| Gowns Double | | | ● | | ● | | ● | | ● |
| Isolation Sign Holder | | ● | ● | ● | ● | | | ● | ● |
| Gloves Double | | ● | ● | | ● | | | | |
| Gloves Triple | ● | Add GD-011 or GD-111 | | ● | | ● | | | ● |
| Gloves Quad | | Add GD-012 | | | | | | ● | |
| Face Masks Earloop Single | ● | ● | ● | ● | ● | ● | ● | ● | |
| Face Masks Earloop Double | | | | | | | | | ● |
| Face Shield Masks | | ● | | | | ● | ● | ● | |
| Cone Shaped Face Masks | | | | | | | ● | | |
| Storage Add On | MP-015, MP-042, MP-085 | MP-015, MP-042, MP-085 | MP-015, MP-042, MP-085 | | MP-015, MP-042, MP-085 | MP-015, MP-042, MP-085 | MP-015, MP-042, MP-085 | MP-015, MP-042, MP-085 | MP-015, MP-042, MP-085 |
| Storage Built in | | GD-111 | | ● | | | ● | | ● |